



# Easi-Scope Wireless™

Magnifies specimens up to  
43x on a 17" computer monitor



**USER GUIDE**



# CONTENTS

<b>BEFORE USE</b>	<b>4</b>	Button menu	13
Important information	4	System settings menu	14
Care and maintenance	4	Deleting an item	15
Warning	4	<b>INSTALLING THE COMPUTER SOFTWARE</b>	<b>17</b>
Product description	5	Software installation	17
Computer requirements	5	<b>USING THE XPLOVIEW SOFTWARE WITH A COMPUTER</b>	<b>18</b>
<b>PRODUCT AT A GLANCE</b>	<b>6</b>	Connecting the device	18
Package contents	6	Starting the xploview software	18
Product overview	7	Selecting the device in the xploview software	19
Product specification	7	Button menu	20
<b>GETTING STARTED</b>	<b>8</b>	Full screen viewing	21
Charging	8	Image rotation / flip	22
Assembling the product	8	System settings menu	23
Turning on the device	9	Device setup	24
Focusing	9	Timed shot setup	24
Capturing	10	Movie setup	24
Turning off the device	10	Save setting	25
<b>INSTALLING THE SMARTPHONE APP</b>	<b>11</b>	Language setting	25
<b>USING THE XPLOVIEW APP WITH A SMARTPHONE</b>	<b>11</b>	Advanced settings	25
Connecting the device	11	Saved files	27
Open the xploview app	12	Uninstalling the xploview software	27

# BEFORE USE

## Important information

Please read this instruction manual carefully before using this product, and retain this instruction for future reference.

Improvements and changes to this text necessitated by typographical errors, or improvements to the software and / or equipments may be made at any time without notice.

## Care and maintenance

- Avoid vibration, shock and pressure e.g. dropping the microscope.
- Keep the device dry and protect it from water or vapour.
- Do not leave your device in a place with extreme high or low temperature.
- Do not touch the device with a wet hand as it may damage the device, or cause an electric shock to the user.
- Do not use or store the device in dusty, dirty areas as its moving parts may be damaged.
- Do not use harsh chemicals, cleaning solvents or strong detergents to clean the device. Wipe it with a soft cloth slightly dampened in a mild soap-and-water solution.

## Warning

- Do not place lighted wireless microscope upon the eye, doing so may cause permanent eye damage.
- Do not attempt to open or dismantle the wireless microscope.

## Product description

This wireless microscope features a 640 x 480 CMOS sensor and can be paired with up to 2 Wi-Fi enabled devices within a 10m range. It is compatible with most Android, IOS, Windows and Mac OS devices.

Simply connect the device to the wireless microscope's WLAN access point to start viewing, capturing and recording live images.

## Computer requirements

Smartphone (Android and IOS)

**WLAN**

Wi-Fi 802.11 b / g / n

Windows based PC

**Compatible operating systems**

Windows 10 (32 bit or 64 bit)  
Windows 8 (32 bit or 64 bit)  
Windows 7 (32 bit or 64 bit)  
Windows Vista (32 bit or 64 bit)  
Windows XP SP2, SP3

**WLAN**

Wi-Fi 802.11 b / g / n

Mac OS based PC

**Compatible operating systems**

Mac OS X 10.5.6 - Mac OS X 10.11.x


**WLAN**

Wi-Fi 802.11 b / g / n

# PRODUCT AT A GLANCE

## Package contents

Wireless microscope



Installation CD



Adaptor




User manual



USB cable



## Product overview



1. Capture button
2. Focusing ring
3. Power button
4. USB port


5. LED lights
6. Adaptor slot
7. Indication light

## Product specification

<b>Connection type</b>	802.11 b / g / n
<b>Effective magnifications (on a 17" monitor)</b>	43x
<b>Effective viewing area</b>	8 x 6 mm
<b>Illumination</b>	Four LEDs
<b>Sensor</b>	CMOS
<b>Size</b>	∅70 x 85 mm
<b>Weight</b>	170 grams

# GETTING STARTED

## Charging




Charge the wireless microscope before using it.

Connect the wireless microscope to a computer with the USB cable provided.


Status	Indication light colour
Charging	Red
Charged	Green

## Assembling the product




Insert the adaptor into the adaptor slot.






Fasten the adaptor by turning it gently in the clockwise direction.

### Turning on the device




Press and hold the power button for 2 seconds.

### Focusing




Adjust the focus of the image manually by rotating the focusing ring.

## Capturing




Press the image capture button to capture an image.

## Turning off the device



Press and hold the power button for 2 seconds to turn off the device.

## INSTALLING THE SMARTPHONE APP




Search for "xploview" on Play Store or App Store on your mobile device and install the app.


## USING THE XPLOVIEW APP WITH A SMARTPHONE

### Connecting the device


Turn on the wireless microscope (see P.9).



Go to "Settings".



Turn on Wi-Fi connection.



Search for and connect to "xploview\_xxxxxx" under the list of available devices. ("xxxxxx" is a 6 digit code unique to your device.)

The default password of the device is 12345678.

## Open the xploview app



### **⚠ Note**

The performance of the device may be affected in an environment crowded by Wi-Fi signal.

User may experience a delay in real-time image display depending on the speed of the Smartphone processor.

## Button menu

The icons on the button menu:



Enter Image / Video playback mode

Picture

Picture playback

Video

Video playback



Capture on screen image



Enter Video mode



Start video recording



Pause video recording




End video recording




Open System Settings Menu (see system settings menu on page 14)


## System settings menu




## Deleting an item



Android:  
Press and hold an item for 2 seconds.



Tab "OK" to confirm.



IOS:  
Swipe the item to the left.




Tab "Delete" to confirm.





# INSTALLING THE COMPUTER SOFTWARE

## Software installation

### Windows based PC

- Insert the supplied application CD to the CD-ROM of the computer.
- Double click the "xploview wireless v3.2.xx.exe" icon <  > located on the driver CD.
- Follow the xploview setup wizard to install the application software for the wireless microscope.

### Mac OS based PC

- Insert the supplied application CD to the CD-ROM of the computer.
- Double click the "xploview wireless v3.2.xx.dmg" icon <  > located on the driver CD.
- Drag the xploview icon <  > into the Applications folder.

# USING THE XPLOVIEW SOFTWARE WITH A COMPUTER

## Connecting the device


Turn on / connect the WiFi adaptor in the Windows or Mac operating system.

Search for and connect to "xploview\_XXXXXX" under the list of available wireless networks. ("XXXXXX" is a 6 digit code unique to your device.)


The default password of the device is 12345678.

## Starting the xploview software

Windows based PC

The xploview software can be launched by double clicking the xploview icon <  > from the desktop, or from the start menu.

Mac OS based PC

The xploview software can be launched by double clicking the xploview icon <  > from the Applications menu.

## Selecting the device in the xploview software

From the "Device" drop-down menu under the "System Settings" menu, choose "xploview Wireless" and click the "Connect" button and click "Apply".



## Button menu

The icons on the button menu :



Open **System Settings Menu** (see system settings menu on page23).



**Capture** on screen **image**.



To start and to stop **Timed Shot**. Images will be captured at a regular interval (see system settings menu on page 24 to adjust frequency and duration).



To start and to stop **Video Recording**.




**Information** of the application software. This information maybe helpful when updating software.



**Shut down** the application software.

## Full screen viewing



To activate the full screen mode, click the full screen button <  > located on the bottom right corner of the xplview application software window.

To exit full screen mode, either double click on the screen, or press the "Esc" button on the keyboard.



## Image rotation / flip


Click <  Rotate > button to rotate or preview screen.




## System settings menu

The first time the xploview software is started, the default settings will be loaded, you may change these settings manually in the system settings menu.

Windows based PC



Mac OS based PC



## Device setup


If the image captured by the wireless microscope was not displayed by default, you can change this by selecting "xploview wireless" from the "Device" drop-down menu and click "Connect".



The screenshot shows a dialog box titled "Device Setup". It contains two dropdown menus. The first is labeled "Device:" and has "xploview Wireless" selected. The second is labeled "Resolution:" and has "640 X 480" selected.

## Timed shot setup

The frequency and duration of automatic images capture can be adjusted under this option.



The screenshot shows a dialog box titled "Timed Shot Setup". It has two rows of controls. The first row is for "Take:" with a text input field containing "60", followed by "shots per" and three radio buttons: "minute" (checked), "hour", and "day". A "(max. 300)" label is to the right. The second row is for "Duration:" with a text input field containing "60", followed by three radio buttons: "minute(s)" (checked), "hour(s)", and "day(s)". A "(max. 999)" label is to the right.

## Movie setup

The resolution of videos you record can be changed from the "Resolution" menu. You can also set a maximum file size for each video.




The screenshot shows a dialog box titled "Movie Setup". It contains a dropdown menu labeled "Resolution:" with "640X480" selected. Below it is a text input field labeled "Maximum file size:" containing "50", followed by "MB (16MB - 4097MB)".




## Save setting

The default location for captured images or videos can be changed under this option.



## Language setting


The language of the xploview software can be changed under this option.




## Advanced settings

By clicking the "More..." button on the right of the system settings menu, you will be able to manually adjust all of the image settings. Note that the settings available may be different, depending on your operating system.

## Windows based PC




## Mac OS based PC



## Saved files

With the xploview application software opened, you can locate the saved files folder by clicking the "More..." button located on the left of the main software window.



## Uninstalling the xploview software

Windows based PC

Select "Uninstall xploview" from the start menu (Start > All Programs > xploview wireless > Uninstall xploview).

Mac OS based PC

Drag the xploview wireless application icon from the "Applications" folder to "Trash".

## Warranty & Support

This product is provided with a one year warranty for problems found during normal usage. Mistreating the Easi-Scope Wireless or opening the unit will invalidate this warranty.

All out of warranty repairs will be subject to a replacement fee.

## Technical Support

Please visit [www.tts-group.co.uk](http://www.tts-group.co.uk) for the latest product information.

Email [feedback@tts-group.co.uk](mailto:feedback@tts-group.co.uk) for technical support.

TTS Group Ltd.  
Park Lane Business Park  
Kirkby-in-Ashfield  
Nottinghamshire,  
NG17 9GU, UK.  
Freephone: 0800 318686 Freefax: 0800 137525

TTS is proud to be part of **RM** plc

## Easi-Scope Wireless™

TTS Product Code:  
EL00470

**WARNING:** Do not dispose of this product in household waste. Hand it over to a collection point for recycling electronic appliances.

**WARNING!** Not suitable for children under 36 months due to small parts – choking hazard. Made in China, on behalf of the TTS Group Ltd.

